## **Technical Datasheet**

## SuperTaq Polymerase

 Cat No:
 020505
 Pack size:
 50U

 020511
 250U
 250U

 020503
 500U
 1kU

 020531
 5kU

Unit Definition: One unit is the amount of enzyme required to incorporate 10nmoles total nucleotides into acid-insoluble material in 30 min at 74 °C.

Assay conditions:

- 25mM TAPS (pH 9.3 at 25 °C) 50mM KCl 2mM MgCl<sub>2</sub> 0.2mM dNTPs each 0.5mg/ml activated salmon sperm DNA 1mM 2-mercaptoethanol
- Storage buffer: 50mM Tris-HCI (pH 8.0) 1mM EDTA 1mM DTT 50% v/v glycerol and stabilizers

10X PCR buffer: 100mM Tris-HCl (pH 9.0) 15mM MgCl<sub>2</sub> 500mM KCl 1% Triton X-100 0.1% w/v stabiliser

QC Assay: This enzyme has been tested for the absence of endonuclease and exonuclease activities. It is >95% pure by SDS gel electrophoresis analysis. Each batch is assayed for amplification by PCR of the 506bp fragment of the *Thermus thermophilus* RNase H gene.

PCR Optimisation:

- 1. If the PCR reaction is not as effective as expected under the normal buffer conditions, try increasing Mg concentration to 1.75mM with 0.3mM dNTPS each final concentration.
- 2. For medium-sized amplifications (1-2kb fragment), it is possible to obtain better amplification by carrying out the elongation at 68 °C rather than at 72 °C.
- 3. Use the same amount of enzyme but at a higher concentration (ie: dilute in a smaller volume). This is very important when doing more difficult PCR like RAPDS, PCR of tailed DNA, or single cell copy PCR. A better result can normally be obtained by using, for example, 0.5µl of 1U/µl enzyme, rather than 1µl of 0.5U/µl.

Storage: Store at –20°C.

<sup>&</sup>lt;sup>1</sup> A higher concentration of SuperTaq is also available for RAPDS and other difficult PCR. Please call for details.