Aqua Screen® *FastExtract*

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(*) Contains chaotropic salt, which is an irritant. Take appropriate laboratory safety measures and wear gloves when handling. Not compatible with disinfecting agents containing bleach.

(†) Contains sodium azide as a preservative. When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

1.2 Stability and Storage

Kit components are stored at room temperature. By following these recommendations, the kit is stable until the expiration date. Following the initial use, the reagents are to be used within 2 months. Prior to first use, add the above mentioned volumes of ethanol (96-100%) to Wash Buffer 1 and 2. Bring the reagents to room temperature (+18°C to +25°C) prior to use. Resuspend possible precipitate in the Lysis Buffer by heating. Buffer may be heated up to 70°C without loss of quality.

1.3 Supplemental Requirements

vacuum pump

micro centrifuge

filtration system, 47 or 50 mm frit

micropipettes and filtered tips

ethanol (96-100%)

Incubator (37 °C for petri dishes, 56 °C for reaction tubes)

2. Application and Test Principle

Aqua Screen® is an *in vitro* test system for the qualitative or quantitative detection of water pathogens. The flexibility of the test system allows for a variety of water sources to be tested. Aqua Screen® utilizes the standard 0.45 µm membranes used routinely with filtration systems, the lysis of the precipitated microorganisms directly on the filtration membrane, the DNA extraction as well as the purification and elution of the DNA in minimal volumes under conditions which are optimal for the subsequent PCR analysis.

3. Test Protocol

3.1 Preparation of Sample Material

Drinking water, bathing or pool water and waste water released from suspended particles can be used as sample materials. The sample conditioning affects the reliability of the test findings and must be in accordance with the guidelines of the ISO Norm 11731. For additional guidelines pertaining to the water sampling procedure, please refer to ISO Norm 5667-1 to 5667-10 or contact your national environmental agency. The water samples should be kept at room temperature, and tested within two days of collection. Preservation with chemical means is not possible as intact microorganisms are needed for the filtration procedure. Water samples contaminated with suspended particles or fixed volatile contents have to be purified by filtration with a folded paper filter. The samples may not be centrifuged for purification. For the testing procedure, at least 100 ml is minimally required, however a sample volume of 1000 ml is recommended. Free DNA of already lysed microorganisms in the sample passes through the membrane filter and is not detectable in the test system. The test does not differentiate between "viable and culturable" and "viable but not culturable" (VBNC) legionella.

3.2 Filtration of the Water Sample

1. After carefully removing a membrane filter from the packaging moisten the filter for a few seconds using the sample material.
2. Insert the membrane filter into the filtration system and filtrate the amount of water required.

3.3 Lysis

Resuspend possible precipitate in the Lysis Buffer, by heating it up to 70 °C if necessary.

1. Transfer the membrane upside down into the incubation disks.
2. Pipette **2 ml** of prewarmed **Lysis Buffer**.
3. Submerge the membrane in the Lysis Buffer by carefully rocking the dish and incubate at 37°C for 30 min.
4. Rinse the membrane and homogenize the lysis buffer by rocking the dish for 10 seconds.
5. Pipette 500 µl of the lysate into the 1.5 ml-incubation tube and incubate for 15 min at 56 °C.
6. Add **250 µl of ethanol** and vortex 5 sec.

3.4 DNA Elution

Aqua
FX

Temperate the Elution Buffer to 56°C.

1. Take one spin column per sample from the kit and stick it into a collection tube. Mark the sample identification on the lit of the spin column and fill **700 µl of the sample lysate** into the spin column.
2. Centrifuge the system for 1 min at 8,000 rpm (6,000 rcf).
3. Place the spin column into a new collection tube, and discard the collection tube containing the lysate.
4. Put the spin column into a new collection tube and add **500 µl of Wash Buffer 1**.
5. Centrifuge the system for 1 min at 8,000 rpm (6,000 rcf) in order to remove the remaining Wash Buffer 1.
6. Place the spin column into a new collection tube and fill the spin column with **500 µl Wash Buffer 2**.
7. Centrifuge the system for 1 min at 8,000 rpm (6,000 rcf) in order to remove Wash Buffer 2.
8. Take the spin column out of the collection tube and dump the containing Wash Buffer 2.
9. Place the spin column again into the used collection tube.

10. Centrifuge for 30 sec at 13,200 rpm (16,100 rcf) in order to remove the remaining Wash Buffer 2.
11. Place the spin column into a sample storage tube, and discard the collection tube containing the Wash Buffer 2.
12. Pipette **60 µl of pre-warmed elution buffer** (56 °C) into the spin column. Secure the sample storage tube and incubate for 5 min at room temperature.
13. Following the incubation, centrifuge the system for 1 min at 8,000 rpm (6,000 rcf) followed by 20 sec. at 13,200rpm (16,100 rcf).
14. Remove the spin column and use the eluate directly for the PCR procedure. The extract is stable for approximately 1 week at +2 °C to +8 °C and can be kept at -18 °C for long term storage.

4. Calculations

With the use of quantitative real-time PCR and a standard curve generated with the titrated DNA standard (Cat. No. 52-0101), the quantity of legionella particles in the extract and subsequently in the sample volume can be determined. If a 10 µl test volume is used, then the total quantity of genomic equivalents present in the water sample can be calculated:

$$\text{genomes/PCR} \times 25,71 \times \frac{1000}{\text{sample volume in [ml]}} = \text{legionella particles per liter}$$

Example: A 500 ml water sample is used for the analysis and extracted in 60 µl. On the basis of the real-time PCR, a result of 30 DNA copies is determined using a test volume of 10 µl.

$$30 \text{ genomes/PCR} \times 25,71 \times 2 = 1542 \text{ legionella particles per liter}$$

In the 500 ml sample water, there were 771 intact legionella particles (viable culturable, viable but not culturable (VBNC), and dead yet still intact legionella). That results in 154 legionella particles per 100 ml of water.

5. Safety Information

Always wear a suitable lab coat, disposable gloves, and protective goggles. The sample-preparation waste contains lysis buffer and wash buffer 1, which can form highly reactive compounds when combined with bleach. **DO NOT add bleach or acidic solutions directly to the sample-preparation waste.** If liquid containing these buffers is spilt, clean with suitable laboratory detergent and water.

The risk phrases applying to lysis buffer are:

- R 22 Harmful if swallowed.
- R 36/38 Irritating to eyes and skin.
- R 52/53 Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

The risk phrases applying to wash buffer 1 are:


- R 22 Harmful if swallowed.
- R 36/38 Irritating to eyes and skin.

The risk and safety phrases applying are:


- S 13 Keep away from food, drink and animal feed.
- S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- S 36 Wear suitable protective clothing.
- S 46 If swallowed, seek medical advice immediately and show container or label.

Process:

1. Filtration:

2. Lysis

Related Products

MB Taq DNA Polymerase
53-0050/-0100/-0200/-0250 MB Taq DNA Polymerase 50/100/200/250 units

Quantification Standards, 100 µl each, 1x10⁶ genomes/µl
52-0101 *Legionella pneumophila* DNA Standard

DNA Remover™
15-2025 DNA Decontamination Reagent, spray bottle 250 ml
15-2200 DNA Decontamination Reagent, refill bottles 4 x 500 ml

Aqua Screen Fast Extract Kit
32-1010/1050/1200 Extraction Kit 10/50/250 extractions

Diagnostic Kits for Conventional PCR

33-1025	<i>Legionella</i> spp. Detection Kit for conv. PCR	25 tests
33-1100	<i>Legionella</i> spp. Detection Kit for conv. PCR	100 tests
33-1250	<i>Legionella</i> spp. Detection Kit for conv. PCR	250 tests
34-1025	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	25 tests
34-1100	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	100 tests
34-1250	<i>Legionella pneumophila</i> Detection Kit for conv. PCR	250 tests

Diagnostic Kits for Real-Time PCR

33-2025	<i>Legionella</i> spp. Detection Kit type 1	25 tests
33-2100	<i>Legionella</i> spp. Detection Kit type 1	100 tests
33-2250	<i>Legionella</i> spp. Detection Kit type 1	250 tests
33-4025	<i>Legionella</i> spp. Detection Kit type 2	25 tests
33-4100	<i>Legionella</i> spp. Detection Kit type 2	100 tests
33-4250	<i>Legionella</i> spp. Detection Kit type 2	250 tests
34-2025	<i>L. pneumophila</i> Detection Kit type 1	25 tests
34-2100	<i>L. pneumophila</i> Detection Kit type 1	100 tests
34-2250	<i>L. pneumophila</i> Detection Kit type 1	250 tests
34-4025	<i>L. pneumophila</i> Detection Kit type 2	25 tests
34-4100	<i>L. pneumophila</i> Detection Kit type 2	100 tests
34-4250	<i>L. pneumophila</i> Detection Kit type 2	250 tests

Appendix

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Trademarks

Aqua Screen is a registered trademark of Minerva Biolabs.