



INTRODUCTION

CIRRUS™ (Product code CIR-RNA-1000-01) is a ‘ready-to-use’ pre-formulated master mix that is specifically designed for fluorogenic probe based reverse transcription Polymerase Chain Reaction (RT-PCR) applications using the 5’ nuclease signalling principle and RNA/DNA targets. It is also suitable for other dual labelled fluorogenic probe methods based on conformational signalling. **CIRRUS™** has been formulated using a proprietary blend of sugars, stabilisers and macromolecules in a freeze drying process providing instant reagent dissolution using only water. This drying process ensures the entire active components remain stable at ambient temperatures and eliminates the requirement for refrigerated reagent transport and storage. Once opened the **CIRRUS™** cake is resistant to high humidity levels for up to one hour making it ideal for automation.

CIRRUS™ contains the core reaction components required to carry out a RT-PCR, this includes: core PCR reaction buffer, dNTP (inc dUTP) nucleotides, Magnesium ions and a high performance rapidly activated HotStart *Taq* DNA polymerase and a thermostable *MMuLV* RNA dependant DNA polymerase. This enzyme provides superior conversion of RNA to cDNA at elevated temperatures, increasing primer stringency and reducing RNA secondary structure. Supplied as a dry reagent, **CIRRUS™** may be reconstituted using water and the addition of target-specific oligonucleotide primers and probes (not supplied) to provide a complete RT-PCR reaction mix.

Highly flexible, the **CIRRUS™** reagent is readily compatible with automated laboratory dispensing systems and may be reconstituted to suit a variety of end-user test requirements. This includes the addition of other buffer solutes and third party enhancers such as Uracil-N-glycosylase (UNG).

CIRRUS™ may be reconstituted in a volume that allows more template or reagents to be added than conventional 2x frozen mixes. This may be useful when using multiple primer sets in a multiplex format, or where sensitivity is needed (more template). Each vial contains an ~7.5% reagent excess and should be reconstituted to the following volumes:

Master Mix concentration req'd	Diluent Volume
10X	135 µL
5X	270 µL
2X	675 µL*

***Typical application**

Example protocols are provided below using **CIRRUS™** reagents. Given the unique versatility of **CIRRUS™**, many

experimental protocols are possible and those provided are intended as a guideline only.

PROTOCOLS

Opening of **CIRRUS™** Vial: Remove the cap. Gently tap the glass vial to settle contents that may have shifted during transit. Remove neoprene stopper and resuspend as described in the following user protocols and vortex mix briefly.

Protocol 1: 50 x 25 µL PCR reactions using 5µl template per reaction.

The **CIRRUS™** vial contents may be reconstituted to a volume of 675 µL to provide a 2X core reaction mix. This is then used to produce a sub master containing primers and probes. This is sufficient mix for 50 x 25 µL (final) PCR reactions as follows:

Reagent	Volume
CIRRUS™ vial <i>(Product code CIR-RNA-1000-01)</i> Resuspended to <u>2x</u> concentration using sterile water (675 µL). Withdraw 625 µL.	625 µL
Oligonucleotide Primers <i>(not supplied)</i> Forward Primer to a final reaction concentration of 0.1 to 1µM	e.g. 12.5-125 µL 10µM stock solution
Reverse Primer to final reaction concentration of 0.1 to 1 µM	e.g. 12.5-125 µL 10µM stock solution
Oligonucleotide Probe <i>(not supplied)</i> Probe (Dual labelled hydrolysis probe / Molecular Beacon) to a final reaction concentration 0.05 to 0.2µM	e.g. 31.25-125 µL of 2µM stock
Diluent/ Nuclease-free water	To final Volume
TOTAL	1000 µL

The reaction mix is then aliquotted (50) in 20 µL volumes in PCR reaction tubes and mixed with 5µL nucleic acid template ready for thermocycling. No template controls may use water to substitute template RNA or DNA.

Thermal Cycling: The following are general recommended thermal cycler settings for **CIRRUS™** reagents. Actual hold temperatures, times and transition rates may vary according to PCR assay type and instrument used.

Phase	Hold	Temp (°C)	Time (s)	Rate (°C/s)
Reverse Transcription	Hold	40-55	60-120	3-10
Enzyme Activation	Hold	95	60	3-10
Amplification	Denature	95	5-10	3-10
	Anneal & Extend	50-65	5-30	3-10



INSTRUCTIONS FOR USE

Protocol 2: 50 x 25 µL RT-PCR reactions using 12.5µL template per reaction.

The **CIRRUS™** vial contents may be reconstituted to a volume of 270 µL to provide a 5X core reaction mix. This is then used to produce a 2x sub master containing primers and probes which can be added to template in equal volume. This is sufficient mix for 50 x 25 µL (final) PCR reactions as follows:

Reagent	Volume
CIRRUS™ vial (Product code CIR-RNA-1000-01) Resuspended to 5x concentration using sterile water (270 µL). Withdraw 250 µL.	250 µL
Oligonucleotide Primers (not supplied) Forward Primer to a final reaction concentration of 0.1 to 1µM Reverse Primer to final reaction concentration of 0.1 to 1 µM	e.g. 12.5-125 µL 10µM stock solution e.g. 12.5-125 µL 10µM stock solution
Oligonucleotide Probe (not supplied) Probe (Dual labelled hydrolysis probe / Molecular Beacon) to a final reaction concentration 0.05 to 0.2µM	e.g. 31.25-125 µL of 2µM stock
Diluent/ Nuclease-free water	To 625 µL final Volume
TOTAL (for 2X reaction mix)	625 µL

The 2X reaction mix is then aliquotted (50) in 12.5 µL volumes in PCR reaction tubes and supplemented with 12.5µL of RNA or DNA template ready for thermocycling. No template controls may use water to substitute template RNA or DNA.

Thermal Cycling: As per protocol 1

ADDITIONAL INFORMATION

CIRRUS™ reagent is formulated to provide 5mM final reaction Magnesium ions. This is optimum for most 5' nuclease assays. This concentration can be increased to suit the requirements of different assays by substituting diluent with additional MgCl₂ solution (not supplied) in the final master mix.

CIRRUS™ reagent is formulated to provide a base performance. This formulation can be adjusted to suit the requirements of different assays/ PCR instruments by substituting diluent with many additional reaction adjuncts.

STORAGE

CIRRUS™ mastermix is supplied as a dried reagent and should be stored in its original packaging at 15-30°C. Once reconstituted, the **CIRRUS™** mastermix will remain stable for 24 hours if stored at 2-8°C. For longer term storage as a conventional mastermix reagent, **CIRRUS™**

mastermix should be resuspended to 2x concentration with a diluent supplemented with 8-16% (v/v) molecular biology grade glycerol and stored at -20°C. Use within 6 months.

TECHNICAL SPECIFICATION

Specification	Dimension
DNA dependant DNA Polymerase	High Performance Hot Start Taq polymerase (derived from <i>Thermus aquaticus</i>)
RNA dependant DNA Polymerase	High performance thermostable RNase H+ recombinant MMuLV
Nucleotides	dNTP containing proprietary dUTP mix
Buffer	Tris, pH 8.8 5' nuclease assay specific salts & enhancers
Magnesium Chloride	5 mM May be increased by user
Storage	15 to 30 °C
Shelf life	18 Months from manufacture date
Dissolution time	< 1s
Volume (final) upon dissolution	User defined >135µL. (Final reaction volume 1350µL)
BSA	Contains Bovine Serum Albumin of USA origin certified BSE free.

MANUFACTURER DETAILS

Fluorogenics LIMITED
227 Tetricus Science Park
DSTL Porton Down
Salisbury
SP4 0JQ
United Kingdom
W: www.fluorogenics.co.uk
E: info@fluorogenics.co.uk
T: +44 2380 972030
F: +44 1722 638022

DISCLAIMER

For Research Use Only, this product may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Fluorogenics Limited's products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Fluorogenics Ltd.

CIRRUS™ contains GoTaq® Hot Start Polymerase manufactured by Promega Corporation for distribution by Fluorogenics Ltd. Licensed to Promega under U.S. Patent Nos. 5,338,671 and 5,587,287 and their corresponding foreign patents.

The practice of the patented polymerase chain reaction (PCR) process may require a license depending on geographical territory as may other processes such as the 5' nuclease assay or Molecular Beacon. Other Fluorogenic processes may also require a license. No other rights are conveyed expressly, by implication or by Estoppel under any patent claim. **CIRRUS™** is a trademark of Fluorogenics Ltd. All other marks are the property of their respective owners. ©2014 Fluorogenics Ltd.



Information/ specifications contained within this leaflet are subject to change without notice. Please check with www.fluorogenics.co.uk for current information.

Product:	RNA Master Mix	Page 2 of 2
Product Code:	CIR-RNA-1000-01	
Document Code-Revision:	FGL-IFU-CIR-1000-RNA-11	
Date of Issue:	20-May-15	