

Cat.No. 317-150

Ribospin™ Seed/Fruit

TOTAL RNA PURIFICATION HANDBOOK

Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.) : tech@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

Hybrid-Q™, Exprep™, Exfection™, Expin™, Exgene™, GenEx™, DirEx™, Hybrid-R™, RiboEx™, Ribospin™, Ribodear™, Allspin™, RiboSaver™, EzClear™, EzSep™, EzPure™, EzPass™, AmpONE™, AmpMaster™, HyperScript™, ProtinEx™, PAGESTA™ are trademarks of GeneAll Biotechnology co., ltd.

© 2015 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in :

GeneAll® Ribospin™ Seed/Fruit (317-150)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

Brief Protocol



Homogenization

Protocol I

100 mg Seed or Fruit

Protocol II

100 mg Grain or Rhizome
(Starch-enriched sample)



Lyse

500 ul Buffer SL
+ 500 ul Buffer ML
+ 10 ul 2-ME
Centrifugation] 1 min, 10.000 xg

500 ul Buffer SL + 5 ul 2-ME
Centrifugation] 1 min, 10.000 xg
300 ul Supernatant + 300 ul Buffer ML



Filter

600 ul lysate (on EzPure filter)
Centrifugation] 1 min, 10.000 xg



Bind

500 ul Supernatant + 250 ul Absolute ethanol (on column)
Centrifugation] 1 min, 10.000 xg



Digest DNA

500 ul Buffer RBW (on column)
Centrifugation] 0.5 min, 10.000 xg
70 ul DNase reaction mixture (on column)
Incubation] 10 min, RT



Wash

500 ul Buffer RBW (on column)
Centrifugation] 0.5 min, 10.000 xg
500 ul Buffer RNW (on column)
Centrifugation] 0.5 min, 10.000 xg
Additional centrifugation] 1 min, 10.000 xg



Elute

50 ul Nuclease-free water (on column)
Centrifugation] 1 min, 10.000 xg



Downstream application



INDEX

| | |
|--|----|
| Index | 05 |
| Kit Contents | 06 |
| Product Specifications | 06 |
| Product Disclaimer | 08 |
| Quality Control | 08 |
| Storage Conditions | 08 |
| Precautions | 08 |
| Product Description | 09 |
| Protocol I. For seed and fruit | 10 |
| Protocol II. For starch-enriched grain and rhizome | 12 |
| Troubleshooting Guide | 14 |
| Appendix I. | 16 |
| Appendix II. | 17 |
| Ordering Information | 18 |

Kit Contents

Cat. No. 317-150

(50 prep)

| Components | Quantity | Storage |
|--|----------|----------------------------|
| Buffer SL | 30 ml | Room temperature (15~25°C) |
| Buffer ML | 30 ml | |
| Buffer RBW | 60 ml | |
| Buffer RNW | 30 ml | |
| Buffer DRB | 5 ml | |
| Nuclease-free water | 15 ml | |
| GeneAll® Column type F (Blue ring) with collection tube | 50 ea | |
| GeneAll® EzPure™ filter column (Yellow) with collection tube | 50 ea | |
| 1.5 ml microcentrifuge tube | 50 ea | |
| Protocol handbook | 1 ea | |
| DNase I | 120 ul | -20°C |

Materials Not Provided

- Reagent : Absolute ethanol, β -mercaptoethanol (ACS grade or better)
- Disposable material : RNase-free pipette tips, Disposable gloves
- Equipment : Microcentrifuge, Equipment for disrupting sample

Product Specifications

Ribospin™ Seed/Fruit RNA mini

| | |
|--|---------------|
| Type | Spin |
| Maximum amount of starting samples | 100 mg / prep |
| No. of preparation | 50 |
| Preparation time | ~ 30 minutes |
| Maximum loading volume of mini spin column | 750 ul |
| Minimum elution volume | 30 ul |

■ Protocol selecting guide for starting sample

| | The list of sample applied with Protocol I | The list of sample applied with Protocol II |
|-----------------|---|--|
| Seeds | <p><i>Capsella bursapastoris</i> (Shepherd's purse) <i>Ulmus davidiana</i> var. <i>japonica</i> (Elm) <i>Daucus carota</i> (Carrot) <i>Raphanus sativus</i> var. <i>sativus</i> (Radish) <i>Zinnia violacea</i> (Garden zinnia) <i>Prunus armeniaca</i> (Apricot tree) <i>Apium graveolens</i> (Celery) <i>Pastinaca sativa</i> (Parsley) <i>Vitis vinifera</i> (Grape tree) <i>Cucurbita</i> spp. (Pumpkin) etc.</p> | <p><i>Phaseolus vulgaris</i> (Kidney bean) <i>Phaseolus radiatus</i> (Mung beans) <i>Triticum aestivum</i> (Wheat) <i>Zea mays</i> (Corn) <i>Setaria italica</i> (Millet) etc. (other starch-enriched grains)</p> |
| Fruits | <p><i>Fragaria ananassa</i> (Strawberry) <i>Malus domestica</i> (Apple) <i>Solanum lycopersicum</i> (Tomato) <i>Musa sapientum</i> L. (Banana) <i>Mangifera indica</i> (Mango) <i>Pyrus serotina</i> (Pear) <i>Citrus unshiu</i> (Mandarin) etc.</p> | |
| Rhizomes | | <p><i>Ipomoea batatas</i> (Sweet potato) <i>Solanum tuberosum</i> (Potato) <i>Dioscorea opposita</i> (Yam) etc. (other starch-enriched rhizomes)</p> |

Product Disclaimer

GeneAll® Ribospin™ Seed/Fruit RNA mini kit is for research use only, and should not be used for drug, household, or other unintended uses. All due care and attention should be taken in every procedure in this handbook. Please consult the Material Safety Data Sheet (MSDS) for information regarding hazards and safe handling practices.

Quality Control

GeneAll® Ribospin™ Seed/Fruit RNA mini kit is manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. For consistency of product, the quality certification process is carried out from lot to lot thoroughly, and only the qualified is approved to deliver.

Storage Conditions

All components of GeneAll® Ribospin™ Seed/Fruit RNA mini kit, except DNase I, should be stored at ambient temperature (15~25°C). DNase I should be stored at -20°C for conservation of enzyme activity.

Storage at cold ambient temperature may cause precipitation in Buffer ML. If precipitate is seen, heat the buffer at 37°C and agitate it for re-solubilization.

All components are stable for 1 year under these conditions.

Precautions

The buffers included in GeneAll® Ribospin™ Seed/Fruit RNA mini kit contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector and follow standard safety precautions.

Product Description

Ribospin™ Seed/Fruit RNA mini kit is designed for easy and convenient isolation of total RNA from difficult plant tissues such as seeds, fruits, and rhizomes. Especially, this kit can remove effectively large quantities of secondary metabolites including polysaccharides and polyphenolic compounds which can lead to inhibition of downstream application.

Ribospin™ Seed/Fruit RNA mini kit provides two different procedures that are available for application of various plant tissues as follows: Protocol I for seed and fruit, Protocol II for starch-enriched grain and rhizome. For efficient RNA purification, this kit offers optimized lysis system according to the sample type and adopts EzPure™ filter column to eliminate impurities simply from lysate. Moreover, contamination of genomic DNA, that causes interference in RNA analysis, can be excluded by on-column DNase I treatment in these procedures.

The purified RNA is suitable for use in various downstream procedures including cDNA synthesis, RT-PCR, or Northern blotting.

For seed and fruit

Before starting

Thaw DNase I enzyme for use on ice.

Prepare DNase I reaction mixture just before step 9.

(DNase I reaction mixture : Mix 2 ul of DNase I with 70 ul of Buffer DRB)

- 1. Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 100 mg of ground sample into a 1.5 ml microcentrifuge tube (not provided).**

Quick and complete pulverization with liquid nitrogen is essential for good result in preparation. The commonly used technique for disruption is grinding with a mortar and pestle, however other method such as bead-beater or blender can be a good alternative.

- 2. Add 500 ul of Buffer SL, 500 ul of Buffer ML, and 10 ul of β -mercaptoethanol to the sample and vortex vigorously for 15 seconds.**

Buffer ML tends to congeal with starch-enriched samples such as grain and rhizome. It is hard to separate supernatant containing RNA from debris. Therefore, if the lysate solidifies after addition of Buffer ML, use "Protocol II" instead that is special procedure for RNA extraction from starch-enriched sample.

- 3. Incubate the mixture for 3 minutes at room temperature.**
- 4. Centrifuge the lysate at 13,000 rpm ($\geq 10,000$ xg) for 1 minute and transfer 600 ul of the supernatant to an EzPure™ filter column (Yellow).**
- 5. Centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute and transfer 500 ul of the pass-through to a new 1.5 ml microcentrifuge tube (not provided).**

Through this step, large cell debris and most of genomic DNAs are filtered on the EzPure™ filter column and small pellet as debris will be formed at the bottom of the collection tube. Be careful not to disturb the pellet when transferring supernatant.

- 6. Add 250 ul of absolute ethanol to the supernatant and mix it well by inversion.**

Do not centrifuge at this step.

After addition of absolute ethanol, precipitates may be visible in the mixture which not affect RNA purification.

- 7. Apply all of the mixture into a mini spin column (Type F, Blue ring) and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute.**

Transfer all solution including any precipitates on the mini spin column.

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 8. Add 500 ul of Buffer RBW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 9. Apply 70 ul of DNase I reaction mixture onto the center of the mini spin column for gDNA digestion. Incubate for 10 minutes at room temperature.**

To make DNase I reaction mixture, prepare 2 ul of DNase I with 70 ul of Buffer DRB per isolation. DNase I is sensitive to physical damage. Therefore, do not mix vigorously. If you want to DNase I treatment in RNA eluate, skip step 8 and 9 and refer to "Appendix I".

- 10. Add 500 ul of Buffer RBW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 11. Add 500 ul of Buffer RNW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 12. Centrifuge at maximum speed for an additional 1 minute to remove residual wash buffer. Transfer the mini spin column to a new 1.5 ml microcentrifuge tube (provided).**

Residual ethanol may interfere with downstream reaction. Make sure that the membrane of column has to be dried completely.

- 13. Add 50 ul of Nuclease-free water to the center of the membrane in mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute.**

To increase the RNA concentration, reduce the volume of elution to 30 ul.

The purified RNA should be put on ice immediately for accurate analysis or stored at -70°C for long-term storage.

For starch-enriched grain and rhizome

Before starting

Thaw DNase I enzyme for use on ice.

Prepare DNase I reaction mixture just before step 10.

(DNase I reaction mixture : Mix 2 ul of DNase I with 70 ul of Buffer DRB)

- 1. Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 100 mg of ground sample into a 1.5 ml microcentrifuge tube (not provided).**

Quick and complete pulverization with liquid nitrogen is essential for good result in preparation. The commonly used technique for disruption is grinding with a mortar and pestle, however other method such as bead-beater or blender can be a good alternative.

- 2. Add 500 ul of Buffer SL and 5 ul of β -mercaptoethanol to the sample and vortex vigorously for 15 seconds.**
- 3. Incubate the mixture for 3 minutes at room temperature.**
- 4. Centrifuge the lysate at 13,000 rpm ($\geq 10,000$ xg) for 1 minute and transfer 300 ul of the supernatant to a new 1.5 ml microcentrifuge tube (not provided).**
- 5. Add 300 ul of Buffer ML to the supernatant and vortex vigorously for 15 seconds and transfer all of the mixture to an EzPure™ filter column (Yellow).**
- 6. Centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute and transfer 500 ul of the pass-through to a new 1.5 ml microcentrifuge tube (not provided).**

Through this step, large cell debris and most of genomic DNAs are filtered on the EzPure™ filter column and small pellet as debris will be formed at the bottom of the collection tube. Be careful not to disturb the pellet when transferring supernatant.

- 7. Add 250 ul of absolute ethanol to the supernatant and mix it well by inversion.**

Do not centrifuge at this step.

After addition of absolute ethanol, precipitates may be visible in the mixture which not affect RNA purification.

- 8. Apply all of the mixture into a mini spin column (Type F, Blue ring) and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute.**

Transfer all solution including any precipitates on the mini spin column.

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 9. Add 500 μ l of Buffer RBW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 10. Apply 70 μ l of DNase I reaction mixture onto the center of the mini spin column for gDNA digestion. Incubate for 10 minutes at room temperature.**

To make DNase I reaction mixture, prepare 2 μ l of DNase I with 70 μ l of Buffer DRB per isolation. DNase I is sensitive to physical damage. Therefore, do not mix vigorously. If you want to DNase I treatment in RNA eluate, skip step 9 and 10 and refer to "Appendix I".

- 11. Add 500 μ l of Buffer RBW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 12. Add 500 μ l of Buffer RNW to the mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 30 seconds.**

After centrifugation, discard the pass-through and re-insert the mini spin column back into the same collection tube.

- 13. Centrifuge at maximum speed for an additional 1 minute to remove residual wash buffer. Transfer the mini spin column to a new 1.5 ml microcentrifuge tube (provided).**

Residual ethanol may interfere with downstream reaction. Make sure that the membrane of column has to be dried completely.

- 14. Add 50 μ l of Nuclease-free water to the center of the membrane in mini spin column and centrifuge at 13,000 rpm ($\geq 10,000$ xg) for 1 minute.**

To increase the RNA concentration, reduce the volume of elution to 30 μ l.

The purified RNA should be put on ice immediately for accurate analysis or stored at -70°C for long-term storage.

Troubleshooting Guide

| Facts | Possible Causes | Suggestions |
|------------------------|---|--|
| Low or no yield | Incorrect use of lysis buffer | According to the sample type, the process for lysis is different. Starch-enriched samples such as grains and rhizomes have to be processed by "Protocol II" for effective RNA extraction. |
| | Too much starting sample | Using too much sample leads to inefficient lysis followed by poor RNA yield. Reduce the amount of starting material. |
| | Insufficient pulverization | For best result, sample should be disrupted completely using proper method. |
| | Too low RNA mass in sample | Some samples have low RNA contents. To increase the RNA concentration in eluate, reduce the volume of elution to 30 ul. |
| RNA degradation | Incorrect treatment of β-mercaptoethanol during lysis | Ensure that the correct volume of β -mercaptoethanol is used in lysis buffer for RNase elimination. The effective amount of β -mercaptoethanol is 1% of the lysis volume. |
| | Improper storage of extracted RNA | The purified RNA should be stored at -70°C for long-term storage. Do not store at -20°C . If possible, perform downstream application immediately for accurate analysis after RNA extraction. |
| | RNase contamination | To prevent RNA degradation, wear gloves during all procedure and use RNase-free products with sterile and disposable plastic ware. |
| | Too old starting sample | After sufficient pulverization of starting material, store the sample properly at -70°C . If possible, perform the procedure of RNA extraction immediately after disruption of sample to decrease RNA degradation. |

■ Troubleshooting Guide

| Facts | Possible Causes | Suggestions |
|--|--|---|
| Clogging of EzPure™ filter column | Solidification of lysate | According to the sample type, the process for lysis is different. If the lysate solidifies during lysis of protocol I, the sample may contain a lot of carbohydrate and polysaccharide. For effective RNA extraction from starch-enriched samples, apply "Protocol II". |
| | High viscosity of lysate | RNA can be sheared in viscous lysate that causes clogging of column. Increase centrifugal g-force and time to solve clogging if necessary. |
| Clogging of type F mini spin column | Low centrifugal force | Increase g-force (< 10,000 xg) and time (~3min). |
| | Opaque or viscous binding mixture | According to the sample type, the lysate mixed with ethanol becomes opaque or viscous. It does not affect RNA purification. However, if column is clogged because of these problems, increase centrifugal g-force and time until all mixture passes through the membrane of mini spin column. |
| DNA contamination of RNA eluate | High DNA mass in sample | Some plant tissues have high DNA contents. In this case, genomic DNA can be included in RNA eluate. To reduce DNA contamination effectively, refer to the appendix I "DNase I treatment in eluate". |
| | Incorrect treatment of DNase I reaction mixture | For sufficient enzymatic reaction, add DNase I reaction mixture onto the center of the membrane in mini spin column. |

DNase I treatment in eluate

Appendix I describes how to use the DNase I (included in this kit) to eliminate contaminating genomic DNA in RNA eluate. For samples containing high DNA contents, this method is strongly recommended. This procedure is more efficient than on-column DNase I treatment.

Protocol

1. Prepare the mixture as below in a 1.5 ml microcentrifuge tube.
 - 50 ul RNA eluate
 - 5 ul Buffer DRB
 - 1 ul DNase I
2. Incubate the mixture for 10 minutes at room temperature.
3. Add 1 ul of 0.25 M EDTA per 50 ul eluate.
4. Inactivate DNase I enzyme at 75°C for 10 minutes.

** For efficient DNase I treatment and clean-up of eluated RNA, use of Riboclear™ plus (Cat.No 313-150) is suggested.*

■ Related product

| Product | Cat.No | Size | Features and Benefits |
|------------------------|---------|---------|---|
| Riboclear™ plus | 313-150 | 50 prep | <ul style="list-style-type: none"> - Preparation time : ~17 minutes - High recovery rate : ~95 % - Stable and consistent yield - Efficient removal of genomic DNA including DNase I - Concentrated RNA eluate using micro column - Complete removal of salt and enzymes - No need of additional material - No use of organic solvents, no ethanol precipitation |

Electrophoresis method for using formaldehyde-agarose gel (Denaturing gel method)

A denaturing agarose gel is routinely used for the assessment of the quality of extracted RNA. The RNA isolated from samples forms secondary structure via intramolecular base pairing. Therefore, it is very difficult to analyze the result of electrophoresis because of migrating inaccuracy. However, the formaldehyde-agarose gel denatures the secondary structure of RNA, making accurate migration.

To confirm the RNA band after electrophoresis, the gel should be transferred to a UV transilluminator. Mainly, two RNA bands are shown. If they are intact, the RNA bands should be sharp and the intensity of upper band should be about twice compared to that of the lower band.

Preparation of denaturing gel

1. Put 1g agarose in 72 ml water and heat to dissolve thoroughly.
2. Cool to 60°C.
3. Add 10 ml of 10X MOPS buffer, 18 ml of 37% formaldehyde, and 1 ul of 10 mg/ml ethidium bromide (EtBr).
4. Mix well then pour the gel into the gel tray and cool to solidify it.
5. Transfer the solidified gel from tray to tank, and add enough 1X MOPS running buffer to cover the gel.

Preparation of RNA sample

1. Make the mixture.

| | |
|--|---------------------|
| - ? ul RNA (up to 20 ug) | - 4 ul formaldehyde |
| - 2 ul 10X MOPS electrophoresis buffer | - 10 ul formamide |
2. Incubate the mixture for 15 minutes at 65°C.
3. Chill the sample for 5 minutes in ice.
4. Add 2 ul of 10X formaldehyde gel-loading dye to the mixture.
5. Load the mixture in a denaturing gel which is covered with a sufficient 1X MOPS electrophoresis buffer.
6. Run the gel and confirm the RNA band on transilluminator. Occasionally, destaining gel in dH₂O for several hours may be needed to increase the visibility of the RNA band.

■ Composition of buffers

| 10X MOPS buffer | 10X formaldehyde gel-loading dye |
|------------------------|----------------------------------|
| - 0.2 M MOPS | - 50 % glycerol |
| - 20 mM sodium acetate | - 10 mM EDTA |
| - 10 mM EDTA | - 0.25 % (w/v) bromophenol blue |
| - pH to 7.0 with NaOH | - 0.25 % (w/v) xylene cyanol FF |

* Caution

When handling of formaldehyde-agarose gel, always use gloves and eye protector to avoid contact with skin and eyes. Especially, formaldehyde and ethidium bromide (EtBr) should be handled in a fume hood.

Ordering Information

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA

| | | | |
|-------------------|-----|---------|-------------|
| Plasmid Rapidprep | 50 | 100-150 | mini / spin |
| | 200 | 100-102 | |

GeneAll® Expres™ for preparation of plasmid DNA

| | | | |
|------------|-------|---------|---------------|
| Plasmid SV | 50 | 101-150 | spin / vacuum |
| | 200 | 101-102 | |
| | 1,000 | 101-111 | |
| Midi | 26 | 101-226 | spin / vacuum |
| | 50 | 101-250 | |
| | 100 | 101-201 | |

GeneAll® Exfection™ for preparation of highly pure plasmid DNA

| | | | |
|--------------------------------|-----|---------|---------------|
| Plasmid LE (Low Endotoxin) | 50 | 111-150 | spin / vacuum |
| | 200 | 111-102 | |
| Midi | 26 | 111-226 | spin / vacuum |
| | 100 | 111-201 | |
| Plasmid EF (Endotoxin Free) | 20 | 121-220 | spin |
| | 100 | 121-201 | |

GeneAll® Expin™ for purification of fragment DNA

| | | | |
|------------|-----|---------|---------------|
| Gel SV | 50 | 102-150 | spin / vacuum |
| | 200 | 102-102 | |
| PCR SV | 50 | 103-150 | spin / vacuum |
| | 200 | 103-102 | |
| CleanUp SV | 50 | 113-150 | spin / vacuum |
| | 200 | 113-102 | |
| Combo GP | 50 | 112-150 | spin / vacuum |
| | 200 | 112-102 | |

GeneAll® Exgene™ for isolation of total DNA

| | | | |
|-----------------|-----|---------|---------------|
| Tissue SV | 100 | 104-101 | spin / vacuum |
| | 250 | 104-152 | |
| | 26 | 104-226 | spin / vacuum |
| | 100 | 104-201 | |
| MAXI | 10 | 104-310 | spin / vacuum |
| | 26 | 104-326 | |
| | 100 | 109-101 | spin / vacuum |
| | 250 | 109-152 | |
| Tissue plus! SV | 26 | 109-226 | spin / vacuum |
| | 100 | 109-201 | |
| | 10 | 109-310 | spin / vacuum |
| | 26 | 109-326 | |

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® Exgene™ for isolation of total DNA

| | | | | |
|-------------------|-----------|---------|---------------|---------------|
| Blood SV | mini | 100 | 105-101 | spin / vacuum |
| | | 250 | 105-152 | |
| | Midi | 26 | 105-226 | spin / vacuum |
| | | 100 | 105-201 | |
| MAXI | 10 | 105-310 | spin / vacuum | |
| | 26 | 105-326 | | |
| | Cell SV | 100 | 106-101 | spin / vacuum |
| | | 250 | 106-152 | |
| MAXI | 10 | 106-310 | spin / vacuum | |
| | 26 | 106-326 | | |
| | Clinic SV | 100 | 108-101 | spin / vacuum |
| | | 250 | 108-152 | |
| MAXI | 26 | 108-226 | spin / vacuum | |
| | 100 | 108-201 | | |
| | MAXI | 10 | 108-310 | spin / vacuum |
| | | 26 | 108-326 | |
| Genomic DNA micro | 50 | 118-050 | spin | |
| | 100 | 117-101 | | |
| Plant SV | mini | 250 | 117-152 | spin / vacuum |
| | | 26 | 117-226 | |
| | Midi | 100 | 117-201 | spin / vacuum |
| | | 10 | 117-310 | |
| MAXI | 26 | 117-326 | spin / vacuum | |
| | 10 | 117-310 | | |
| Soil DNA mini | mini | 50 | 114-150 | spin |
| Stool DNA mini | mini | 50 | 115-150 | spin |
| Viral DNA / RNA | mini | 50 | 128-150 | spin |

GeneAll® GenEx™ for isolation of total DNA

| | | | | |
|---------------|-----|---------|----------|----------|
| GenEx™ Blood | Sx | 100 | 220-101 | solution |
| | | 500 | 220-105 | |
| Lx | 100 | 220-301 | solution | |
| | 100 | 221-101 | | |
| GenEx™ Cell | Sx | 500 | 221-105 | solution |
| | | 100 | 221-301 | |
| GenEx™ Tissue | Sx | 100 | 222-101 | solution |
| | | 500 | 222-105 | |
| | Lx | 100 | 222-301 | solution |
| | | 100 | 222-301 | |

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® GenEx™ for isolation of total DNA

| | | | | |
|--------------------|----|-----|---------|----------|
| GenEx™ Plant | Sx | 100 | 227-101 | solution |
| | Mx | 100 | 227-201 | |
| | Lx | 100 | 227-301 | |
| GenEx™ Plant plus! | Sx | 100 | 228-101 | solution |
| | Mx | 50 | 228-250 | |
| | Lx | 20 | 228-320 | |

GeneAll® DirEx™ series

for preparation of PCR-template without extraction

| | | | | |
|---------------------------|--|------|---------|----------|
| DirEx™ | | 100 | 250-101 | solution |
| DirEx™ Fast-Tissue | | 96 T | 260-011 | solution |
| DirEx™ Fast-Cultured cell | | 96 T | 260-021 | solution |
| DirEx™ Fast-Whole blood | | 96 T | 260-031 | solution |
| DirEx™ Fast-Blood stain | | 96 T | 260-041 | solution |
| DirEx™ Fast-Hair | | 96 T | 260-051 | solution |
| DirEx™ Fast-Buccal swab | | 96 T | 260-061 | solution |
| DirEx™ Fast-Cigarette | | 96 T | 260-071 | solution |

GeneAll® RNA series

| | | | | |
|------------------------|------|-----|---------|----------|
| RiboEx™ | mini | 100 | 301-001 | solution |
| | | 200 | 301-002 | |
| Hybrid-R™ | mini | 100 | 305-101 | spin |
| Hybrid-R™ Blood RNA | mini | 50 | 315-150 | spin |
| Hybrid-R™ miRNA | mini | 50 | 325-150 | spin |
| RiboEx™ LS | mini | 100 | 302-001 | solution |
| | | 200 | 302-002 | |
| Riboclear™ | mini | 50 | 303-150 | spin |
| Riboclear™ plus! | mini | 50 | 313-150 | spin |
| Ribospin™ | mini | 50 | 304-150 | spin |
| Ribospin™ vRD | mini | 50 | 302-150 | spin |
| Ribospin™ vRD plus! | mini | 50 | 312-150 | spin |
| Ribospin™ vRD II | mini | 50 | 322-150 | spin |
| Ribospin™ Plant | mini | 50 | 307-150 | spin |
| Ribospin™ Seed / Fruit | mini | 50 | 317-150 | spin |
| Allspin™ | mini | 50 | 306-150 | spin |
| RiboSaver™ | mini | 100 | 351-001 | solution |

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® AmpONE™ for PCR amplification

| | | | | |
|---|----------|-------------------|---------|-------------|
| Taq DNA polymerase | | 250 U | 501-025 | (2.5 U/μL) |
| | | 500 U | 501-050 | |
| | | 1,000 U | 501-100 | |
| α-Taq DNA polymerase | | 250 U | 502-025 | (2.5 U/μL) |
| | | 500 U | 502-050 | |
| | | 1,000 U | 502-100 | |
| α-Pfu DNA polymerase | | 250 U | 504-025 | (2.5 U/μL) |
| | | 500 U | 504-050 | |
| | | 1,000 U | 504-100 | |
| Fast-Pfu DNA polymerase | | 250 U | 505-025 | (2.5 U/μL) |
| | | 500 U | 505-050 | |
| | | 1,000 U | 505-100 | |
| Hotstart Taq DNA polymerase | | 250 U | 531-025 | (2.5 U/μL) |
| | | 500 U | 531-050 | |
| | | 1,000 U | 531-100 | |
| Taq Premix | 96 tubes | 20 μL | 521-200 | lyophilized |
| | | 50 μL | 521-500 | |
| | | 20 μL | 526-200 | |
| 50 μL | 526-500 | | | |
| α-Taq Premix | 96 tubes | 20 μL | 522-200 | lyophilized |
| | | 50 μL | 522-500 | |
| | | 20 μL | 527-200 | |
| 50 μL | 527-500 | | | |
| HS-Taq Premix | 96 tubes | 20 μL | 525-200 | solution |
| | | 50 μL | 525-500 | |
| | | 20 μL | 520-200 | |
| α-Pfu Premix | 96 tubes | 50 μL | 523-500 | solution |
| Taq Premix (w/o dye) | 96 tubes | 20 μL | 524-200 | lyophilized |
| dNTPs mix | | 500 μL | 509-020 | 2.5 mM each |
| dNTPs set (set of dATP, dCTP, dGTP and dTTP) | | 1 ml x 4 tubes | 509-040 | 100 mM |

* Each dNTPs is available

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® AmpMaster™ for PCR amplification

| | | | |
|-------------------|-------------------|---------|----------|
| Taq Master mix | 0.5 ml x 2 tubes | 541-010 | solution |
| | 0.5 ml x 10 tubes | 541-050 | solution |
| α-Taq Master mix | 0.5 ml x 2 tubes | 542-010 | solution |
| | 0.5 ml x 10 tubes | 542-050 | solution |
| HS-Taq Master mix | 0.5 ml x 2 tubes | 545-010 | solution |
| | 0.5 ml x 10 tubes | 545-050 | solution |
| α-Pfu Master mix | 0.5 ml x 2 tubes | 543-010 | solution |
| | 0.5 ml x 10 tubes | 543-050 | solution |

GeneAll® HyperScript™ for Reverse Transcription

| | | | |
|---|------------------|---------|----------|
| Reverse Transcriptase | 10,000 U | 601-100 | solution |
| RT Master mix | 0.5 ml x 2 tubes | 601-710 | solution |
| RT Master mix with oligo (dT) ₂₀ | 0.5 ml x 2 tubes | 601-730 | solution |
| RT Master mix with random hexamer | 0.5 ml x 2 tubes | 601-740 | solution |
| RT Premix | 96 tubes, 20 μl | 601-602 | solution |
| RT Premix with oligo (dT) ₂₀ | 96 tubes, 20 μl | 601-632 | solution |
| RT Premix with random hexamer | 96 tubes, 20 μl | 601-642 | solution |
| One-step RT-PCR Master mix | 0.5 ml x 2 tubes | 602-110 | solution |
| One-step RT-PCR Premix | 96 tubes, 20 μl | 602-102 | solution |
| First strand Synthesis Kit | 50 reaction | 605-005 | solution |
| ZymAll™ RNase Inhibitor | 10,000 U | 605-010 | solution |
| ZymAll™ RNase Inhibitor | 4,000 U | 605-004 | solution |

| Products | Scale | Size | Cat. No. | Type |
|----------|-------|------|----------|------|
|----------|-------|------|----------|------|

GeneAll® Protein series

| | | | |
|--|-----------------|---------|----------|
| ProtinEx™ Animal cell / tissue | 100 ml | 701-001 | solution |
| PAGESTA™ Reducing 5X SDS-PAGE Sample Buffer | 1 ml x 10 tubes | 751-001 | solution |

NOTE

NOTE

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Mailing Address

GeneAll Biotechnology Co., LTD.
GeneAll Bldg., 303-7 Dongnam-ro,
Songpa-gu, Seoul, Korea 138-859

Ordering information

Tel: 82-2-407-0096

Fax: 82-2-407-0779

E-mail: sales@geneall.com

Technical information

Tel: 82-2-407-0096

Fax: 82-2-407-0779

E-mail: tech@geneall.com

Customer & Technical Support

Do not hesitate to ask us any question.
We thank you for any comment or advice.



www.geneall.com

GeneAll Bldg., 303-7 Dongnam-ro
Songpa-gu, Seoul, Korea 138-859

E-mail : sales@geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

©2015 GeneAll Biotechnology, All right reserved

Edited by BR
Designed by Park So Yoon