

QuickExtract – Rapid and efficient extraction of PCR-ready genomic DNA from plant and seed samples

Alex Kolb, Luke Linz, Hank Daum, Jonathan Curry, Dietrich Köster and Julie Kramer

LGC, Biosearch Technologies: Unit 1-2 Trident Industrial Estate, Pindar Road, Hoddesdon, EN11 0WZ, UK

Corresponding author: Dietrich.Koester@lgcgroup.com

Introduction

The screening of large populations and the associated challenges in automation, chemistry and cost are major scientific and business drivers in plant breeding. The QuickExtract™ Plant DNA Extraction

Solution can be used to rapidly and efficiently extract PCR-ready genomic DNA from most plant leaf samples using a simple, one-tube protocol that takes only 8 minutes.

This study shows the successful use of QuickExtract. The influence of grinding on the quality of genotyping results was investigated for multiple plant species. Leaf and seed material were used for testing.

PCR based genotyping methods such as KASP® or BHQplus® Probes were tested. In addition, it is shown that the QuickExtract lysate can be stored and re-used even days later.

Materials and methods

A. Determination of the impact of grinding

Whole seeds from wheat, tomato, hulled sunflower and pepper were either ground or incubated whole with QuickExtract. 100 µL of QuickExtract was added to tomato and pepper seeds, 200 µL was added to wheat seeds and 300 µL was added to sunflower seeds. Samples were processed according to the QuickExtract workflow in Fig 1, and diluted 4 and 16 fold prior to PCR amplification with KASP on the IntelliQube®.

B. Evaluate BHQplus and KASP chemistries for relevant crop types

QuickExtract Solution was added to 24 samples of seed and leaf material of multiple crop types according to Fig 1 and Table 1. Samples were incubated at 65 °C for 8 mins, diluted 2-8 fold (Fig 2) and PCR amplified on the IntelliQube with KASP and BHQ chemistries (see Fig 3).




Figure 1: Overview about the QuickExtract workflow.

After cell lysis, inhibitory cell components are degraded, whereas the released DNA can be inserted directly into the PCR after further heat incubation.

Crop	Tissue Type	Processing	Grinding required	Mass (mg)	QuickExtract Buffer added (µL)
Corn	Seed	Chipped	Optional	~20 mg	200
	Leaf	Punches	Optional	~100 mg	400
Wheat	Seed	Whole	Required	~30 mg	200
	Leaf	Out	Required	~5-10 mg	100
Rapeseed	Seed	Whole	Optional (Grinding slightly better)	~45 mg	400
	Leaf	Punches	Grinding not tested	~5-10 mg	100
Soy	Seed	Half	Grinding not tested	~80 mg	200
	Leaf	Whole	Optional	~3 mg	100
Tomato	Seed	Punches	Optional	5 punches	200
	Leaf	Punches	Optional	15 punches	400
Pepper	Seed	Whole	No Grinding	~8 mg	100
	Leaf	Punches	No Grinding	5 punches	200
Cotton	Seed	Whole	Required	~100 mg	200
	Leaf	Hulled	Required	~50 mg	200
Sunflower	Seed	Whole	Not tested	2 punches	200
	Leaf	Punches	Not tested	2 punches	200

Table 1: Sample and pre-treatment of validated crops, seeds or leaves, grinding requirement conditions, starting weight, and volume of QuickExtract required.

C. Stability of extracted DNA for at least 4 weeks at 4 °C

Corn and tomato leaf samples were extracted with QuickExtract. The DNA was PCR amplified using 2 different KASP assays for each crop type. The sample plate was stored at 4 °C and the experiment was repeated 1 month later (see Fig 4).

Results

A. Determination of the impact of grinding




Figure 2

B. Evaluate BHQplus and KASP chemistries for relevant crop types




Figure 3

C. Stability of extracted DNA for at least 4 weeks at 4 °C




Figure 4

Conclusion

The results show that grinding needs to be determined for each sample used. QuickExtract allows genotyping of tomato and sunflower seeds without grinding. Grinding was required for amplification with wheat seeds, whereas grinding was inhibitory to PCR for pepper seeds. Similar results were obtained with leaf material from these crops (data not shown). Efficient cleaning of the seeds due to coating with chemicals need to be considered.

QuickExtract can be used to extract DNA from leaf and seed material from multiple types of plants for endpoint PCR applications (KASP and BHQplus probes). The results show whether to grind the sample and sample mass to QuickExtract volume must be empirically determined for each crop and sample type. Volume of the processed sample to use for amplification must also be optimised.

Good genotyping results were obtained with the QuickExtracted samples even after storage for 4 weeks at 4 °C.

Summary

- PCR-quality DNA may be extracted from leaf or seed material from multiple plant species.
- The processed sample is stable if stored at 4 °C for 1 month or for several years when stored at -20 °C.
- Prior to large scale processing of multiple samples the following should be considered:
 - Grinding the leaf or seed material prior to processing
 - Optimising the ratio of sample mass to QuickExtract Solution volume
 - Diluting the processed sample before downstream processing

Figure 2: The effect of grinding on seed material.

Wheat, tomato, sunflower and pepper seeds were either ground (red circles) or extracted whole (blue circles) with QuickExtract. Samples were diluted 4 fold (top) or 16 fold (bottom) prior to PCR amplification with KASP on the IntelliQube. Water was used as negative controls (NTC).

Figure 3: Evaluate BHQplus and KASP chemistries for relevant crop types Cluster plots for 7 crops (leaf or seed samples), using KASP and BHQplus Probes The DNA in the QuickExtract Plant DNA Extraction Solution may be further purified using either beadex® or the MasterPure™ DNA purification kit (data not known).

Figure 4: Stability of DNA stored at 4 °C. DNA extracted from corn and tomato with QuickExtract was PCR amplified (original) with two KASP assays and then stored at 4 °C for 1 month. The samples were then again PCR amplified against the same two assays (1 month).

Lucigen®

BIOSEARCH
TECHNOLOGIES

GENOMIC ANALYSIS BY LGC

Integrated tools.
Accelerated science.

Twitter Facebook LinkedIn @LGCBiosearch lucigen.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or by any information storage and retrieval system, without the written permission of the copyright holder. © LGC Limited, 2016. All rights reserved. GEN0000000101 P-16102017

For research or investigational use only