

# GeneReleaser®

## DNA/RNA Releasing Agent

### **GeneReleaser® Protocol for Cerebrospinal Fluid and Synovial Fluid**

GeneReleaser® is composed of proprietary polymeric materials that quickly facilitate the release of nucleic acids from cells in a form suitable for PCR. By segregating inhibitors that are released during lysis as well as preservation agents that may interfere with amplification, GeneReleaser® provides amplifiable nucleic acids from minute amounts of material. The protocols for CSF and synovial fluid are developed by GeneReleaser® users and have not been validated by BioVentures.

### **GeneReleaser® Protocol for Cerebrospinal Fluid (microbial organisms)**

The following are made as guidelines and may require modifications to accommodate your application. Spiking non-infected CS fluid with a model organism may be required to establish suitable protocols. The following protocol works for microbial organisms and intracellular organisms.

1. Dilute 100µl CSF with 200µl PBS in a standard amplification tube.
2. Centrifuge for 5 minutes to pellet cells and bacteria.
3. Discard the supernatant.
4. Wash the pellet with 100µl of 1X PCR buffer by vortexing, centrifuging 1 minute and discarding the supernatant.
5. Thoroughly resuspend the contents of the GeneReleaser® tube by inverting 10-20 times or vortexing briefly.
5. Add 20µl of GeneReleaser® to cell pellet (usually around 10µl pellet volume) and vortex vigorously to resuspend the cell pellet.
4. Place samples onto thermal cycler, with a heated lid, with the following GeneReleaser® program:

Step	Temperature	Time
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.
8.	65°C	60 sec.
9.	80°C	hold

5. Once program is completed, sample is ready to use as PCR template.
6. Add appropriate volume of mastermix (80µl of a 1.25X master mix containing all components for the amplification)
7. Perform amplification reaction according to your optimized protocol.

\* GeneReleaser® treatment can also be performed in the microwave, see procedure on page 2.

### **GeneReleaser® Protocol for Cerebrospinal Fluid (extracellular virus)**

The following are made as guidelines and may require modifications to accommodate your application. Spiking non-infected CS fluid with a model organism may be required to establish suitable protocols.

1. Place 1µl of CSF into the bottom of a 0.2ml-0.5ml standard amplification tube.
2. Resuspend the GeneReleaser® mixture by vortexing 2-3 seconds or inverting 5-10 times.
3. Add 20µl of GeneReleaser® suspension to the 1µl of CSF in the PCR tube and tightly close the tube lid. Pulse vortex briefly to mix.
4. Place samples onto thermal cycler, with a heated lid, with the following GeneReleaser® program:

Step	Temperature	Time
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.

8.	65°C	60 sec.
9.	80°C	hold

- Once program is completed, sample is ready to use as PCR template.
- Add appropriate volume of mastermix (80µl of a 1.25X master mix containing all components for the amplification)
- Perform amplification reaction according to your optimized protocol.

\* GeneReleaser® treatment can also be performed in the microwave, see procedure on page 2.

#### **GeneReleaser® Protocol for Synovial Fluid**

The following are made as guidelines and may require modifications to accommodate your application. Spiking non-infected synovial fluid with a model organism may be required to establish suitable protocols.

- Dilute 100µl synovial fluid with 200µl PBS in a standard amplification tube.
- Centrifuge for 5 minutes to pellet cells and bacteria.
- Discard the supernatant.
- Wash the pellet with 100µl of 1X PCR buffer by vortexing, centrifuging 1 minute and discarding the supernatant.
- Thoroughly resuspend the contents of the GeneReleaser® tube by inverting 10-20 times or vortexing briefly.
- Add 20µl of GeneReleaser® to cell pellet (usually around 10µl pellet volume) and vortex vigorously to resuspend the cell pellet.

- Place samples onto thermal cycler, with a heated lid, with the following GeneReleaser® program:

Step	Temperature	Time
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.
8.	65°C	60 sec.
9.	80°C	hold

- Once program is completed, sample is ready to use as PCR template.
- Add appropriate volume of mastermix (80µl of a 1.25X master mix containing all components for the amplification)
- Perform amplification reaction according to your optimized protocol.

\* GeneReleaser® treatment can also be performed in the microwave, see procedure below.

#### **Microwave Lysis Protocol:**

We have found that the microwave treatment of specimens affords a rapid sample preparation and facilitates the amplification of the more intractable types of specimens.

##### **A. Evaluation of microwave**

Perform the following experiment to determine the optimal conditions for your tubes and microwave.

- Place 40µl DI water in the same type of tube that you will be using for GeneReleaser® treatment.
- Overlay each tube with mineral oil to prevent evaporation.
- Close the tubes, place in microwave safe rack (polyethylene or propylene) and heat on high for 5 minutes.
- If any caps pop or tubes distort in any manner, then place a separate beaker in the microwave with 150ml of room temperature DI water and repeat the above 3 steps, the beaker of water serves as a heat ballast.
- If tubes open or distort, reduce the power by 10% increments and increase time by 1-minute increments repeating step 4 until tubes no longer open or distort.

**Note: Make sure the racks used in this procedure are MICROWAVE SAFE! The 0.2ml tube racks provided by Perkin Elmer for use with their 9600 and 2400 instruments are NOT microwave compatible, they will melt.**

##### **B. Microwave Protocol**

- Perform microwave procedure above for time and power conditions

2. Place 1µl of specimen with 20µl of GeneReleaser® into either a 0.5ml PCR tube or 1.5ml tube.
  3. Vortex the tubes containing specimen and GeneReleaser® for ~10 seconds.
  4. Overlay with mineral oil to prevent samples from evaporating.
  5. Place the closed tubes in a microwave safe polyethylene or propylene rack. Make sure that the lids are loosely closed. If lids are closed too tightly tubes could rupture.
  6. Place the rack in a microwave oven and heat at maximum power setting (setting should be based on the microwave evaluation results) for 5-7 minutes. Typically, 5 minutes if wattage is 900 or higher and 7 minutes if wattage is 500.
  7. Remove rack from microwave and centrifuge the tubes at 5000xg for 5 minutes. After centrifuging samples, remove supernatant and use as DNA template.
  8. Perform the amplification reaction.
- 

#### References:

B E Ley,<sup>a</sup> C J Linton,<sup>b</sup> S Longhurst,<sup>a</sup> H Jalal,<sup>b</sup> M R Millar. Eubacterial approach to the diagnosis of bacterial infection. *Archives of Disease in Childhood* 1997;77:148-149.

Mattu R, Sorbara L, Filie AC, Little R, Wilson W, Raffeld M, Abati A. Utilization of polymerase chain reaction on archival cytologic material: a comparison with fresh material with special emphasis on cerebrospinal fluids. *Modern Pathology* 2004;17: 1295–1301.

If you have further questions, please contact us at [support@bioventures.com](mailto:support@bioventures.com). You may also contact us via fax (1-877-286-0330) or through our toll-free phone number (1-877-852-7846).