

GeneReleaser®

DNA/RNA Releasing Agent

GeneReleaser® Protocol for Mouse Tail

GeneReleaser® is composed of proprietary polymeric materials that quickly facilitate the release of nucleic acids from cells in a form suitable for PCR. By segregating inhibitors that are released during lysis and preservation agents that may interfere with amplification, GeneReleaser® provides amplifiable nucleic acids from minute amounts of material. This manual offers 2 alternate protocols, one tissue homogenization and one Proteinase K digestion.

GeneReleaser® Protocol for Mouse Tail with Tissue Homogenization

Part I: Mouse Tail Homogenization

1. Cut a 1mm thick section from the either fresh or frozen mouse tail that has been rinsed with sterile water to remove any surface contamination and place the section into the bottom of a 1.5 ml tube.
2. Add 25µl of 1X TE to the tube containing the sectioned tail tissue.
3. Mince the section of tissue by thrusting a pestle against the tissue and twisting the pestle to compress the tissue against the walls of the tube. Ten thrusts with the pestle are sufficient.
4. Place the tube containing the homogenate at 4°C as each specimen is homogenized.

Part II: GeneReleaser® Treatment

5. Transfer 1µl of the tissue homogenate obtained above into a 0.5 ml standard amplification tube. NOTE: flick the tube containing the homogenate 3 times to mix before transfer.
6. Resuspend the GeneReleaser® mixture by vortexing 2-3 seconds or inverting 5-10 times.
7. Add 20ul of GeneReleaser® suspension to the 1ul of homogenate in the PCR tube and tightly close the tube lid.
8. Pulse vortex briefly to mix.
9. Place samples onto thermal cycler, with a heated lid, with the following GeneReleaser® lysis program:

Step	Temperature	Time
1.	65°C	30 sec.
2.	8°C	30 sec.
3.	65°C	90 sec.
4.	97°C	180 sec.
5.	8°C	60 sec.
6.	65°C	180 sec.
7.	97°C	60 sec.
8.	65°C	60 sec.
9.	80°C	hold

10. Once program is completed, centrifuge sample tubes at 5,000 x g for 1 minute.
11. Remove supernatant to use as template for PCR. Recommend using 1-10ul of supernatant per 20-100ul amplification reaction. Alternatively, PCR can be performed directly in the GeneReleaser® treatment tube; add amplification reagents for a final volume of 100ul.
12. Perform amplification reaction according to your optimized protocol.

* GeneReleaser® treatment can alternatively be performed in a microwave, see page 2 for the protocol.

GeneReleaser® Protocol for Mouse Tail with Proteinase K Digestion

1. Cut a 1mm thick slice of mouse tail which has been washed with sterile water to remove surface contamination.
2. Place this section of tail in a 0.5ml PCR tube
3. Add 50µl 1XTE.
4. Resuspend the GeneReleaser® mixture by vortexing 2-3 seconds or inverting 5-10 times.
5. Add 20µl GeneReleaser®.
6. Add 2µl Proteinase K (conc.14-15 mg/ml) and mix well.
7. Digest on thermal cycler 1-3 hours @ 55°C.
8. Vortex to resuspend.
9. Heat inactivate 95°C for 10 minutes. (This is critical - DO NOT SHORTEN TIME OR LOWER TEMPERATURE - IF PK IS NOT COMPLETELY INACTIVATED IT WILL DIGEST TAQ.)
10. Centrifuge 5 minutes @ 10,000xg.

11. Carefully transfer supernatant to fresh 1.5ml sterile screw cap tube to use as template for PCR.
12. Use 1 μ l, 2.5 μ l, and 5 μ l of supernatant as template for 3-100 μ l reactions. This is a range finding step. Thereafter use whichever performed best.
13. Perform amplification reaction according to your optimized protocol.

Microwave Lysis Protocol:

We have found that the microwave treatment of specimens affords a rapid sample preparation and facilitates the amplification of the more intractable types of specimens.

A. Evaluation of microwave

Perform the following experiment to determine the optimal conditions for your tubes and microwave.

1. Place 40 μ l DI water in the same type of tube that you will be using for GeneReleaser[®] treatment.
2. Overlay each tube with mineral oil to prevent evaporation.
3. Close the tubes, place in microwave safe rack (polyethylene or propylene) and heat on high for 5 minutes.
4. If any caps pop or tubes distort in any manner, then place a separate beaker in the microwave with 150ml of room temperature DI water and repeat the above 3 steps, the beaker of water serves as a heat ballast.
5. If tubes open or distort, reduce the power by 10% increments and increase time by 1-minute increments repeating step 4 until tubes no longer open or distort.

Note: Make sure the racks used in this procedure are MICROWAVE SAFE! The 0.2ml tube racks provided by Perkin Elmer for use with their 9600 and 2400 instruments are NOT microwave compatible, they will melt.

B. Microwave Protocol

1. Perform microwave procedure above for time and power conditions
2. Place 1 μ l of specimen with 20 μ l of GeneReleaser[®] into either a 0.5ml PCR tube or 1.5ml tube.
3. Vortex the tubes containing specimen and GeneReleaser[®] for ~10 seconds.
4. Overlay with mineral oil to prevent samples from evaporating.
5. Place the closed tubes in a microwave safe polyethylene or propylene rack. Make sure that the lids are loosely closed. If lids are closed too tightly tubes could rupture.
6. Place the rack in a microwave oven and heat at maximum power setting (setting should be based on the microwave evaluation results) for 5-7 minutes. Typically, 5 minutes if wattage is 900 or higher and 7 minutes if wattage is 500.
7. Remove rack from microwave and centrifuge the tubes at 5000xg for 5 minutes. After centrifuging samples, remove supernatant and use as DNA template.
8. Perform the amplification reaction

References:

Kamsler A, Daily D, Hochman A, Stern N, Shiloh Y, Rotman G, Barzilai A. Increased Oxidative Stress in Ataxia Telangiectasia Evidenced by Alterations in Redox State of Brains from Atm-deficient Mice. *Cancer Research* 2001 March; 61:1849–1854.

Oravec-Wilson K, Philips ST, Yilmaz OH, Ames HM, Li L, Crawford BD, Gauvin AM, Lucas PC, Sitwala K, Downing JR, Morrison SJ, Ross TS. Persistence of Leukemia-Initiating Cells in a Conditional Knockin Model of an Imatinib-Responsive Myeloproliferative Disorder. *Cancer Cell* 2009 August 4; 16(2): 137–148.

Sabo T, Lomnitski L, Nyska A, Ben S, Maronpot R, Shoham E, Roses A, Michaelson DM. Susceptibility of transgenic mice expressing human apolipoprotein E to closed head injury: The allele E3 is neuroprotective whereas E4 increases fatalities. *Neuroscience* 2000 November; 4(30):879-884

Stern N, Hochman A, Zemach N, Weizman N, Hammel I, Shiloh Y, Rotman G, Barzilai A. Accumulation of DNA Damage and Reduced Levels of Nicotine Adenine Dinucleotide in the Brains of Atm-deficient Mice. *The Journal of Biological Chemistry* 2002 January; 277:602-608.

If you have further questions, please contact us at support@bioventures.com. You may also contact us via fax (1-877-286-0330) or through our toll-free phone number (1-877-852-7846).