**ZellShield™****The Cell Culture Contamination Preventive****Background**

Zell Shield™ is ready-to-use for protecting cell cultures from a broad range of common contaminants. The extremely powerful combination of innovative antibiotics is blocking the bacterial DNA synthesis, the protein biosynthesis or shows macrolidic properties. These bactericidal features keep contaminants permanently away from the culture. The use of old fashioned antibiotics like penicillin/streptomycin, nystatin or amphotericin B is not longer required! Instead, ZellShield™ can be used for routine addition to cell culture media to prevent contamination with mycoplasma, Gram+ or Gram- bacteria and fungi including yeasts.

Precaution and Storage Information

For research use only! Avoid contact with skin. Don't swallow. ZellShield™ is shipped on cool packs. The product should be stored at 4 °C for immediate use and below -18 °C for long term storage. The product needs to be protected from light.

Instructions for Use

The easy applicable 100x aqueous solution is sterile filtrated and ready-to-go. 1 ml is sufficient for 100 ml of culture medium. One 50 ml bottle will treat 5 liter of culture medium.

ZellShield™ contains a fluoroquinolone and can be used in combination with other antibiotics required as marker. Other components of the product are not applied as genetic markers. The additional use of Pen/Strep is not required or recommended but generally possible. ZellShield™ can also be added directly to any size animal cell culture. Remove and replace fresh ZellShield™ with every time the culture medium is changed but after 7 days of incubation at 37 °C to achieve maximum protection. ZellShield™ can be added to freshly passaged and trypsinated cells or to already precultivated cells by medium change or direct application. The reagent should be refreshed by medium change.

Catalog number: 13-0050, 13-0150

release 08.09

Minerva Biolabs GmbH · Landhausring 7 · 12683 Berlin · Germany
Tel.: +49-30-20004370 · FAX: +49-30-20004379 · E-mail: info@minerva-biolabs.com

**ZellShield™****The Cell Culture Contamination Preventive****Background**

Zell Shield™ is ready-to-use for protecting cell cultures from a broad range of common contaminants. The extremely powerful combination of innovative antibiotics is blocking the bacterial DNA synthesis, the protein biosynthesis or shows macrolidic properties. These bactericidal features keep contaminants permanently away from the culture. The use of old fashioned antibiotics like penicillin/streptomycin, nystatin or amphotericin B is not longer required! Instead, ZellShield™ can be used for routine addition to cell culture media to prevent contamination with mycoplasma, Gram+ or Gram- bacteria and fungi including yeasts.

Precaution and Storage Information

For research use only! Avoid contact with skin. Don't swallow. ZellShield™ is shipped on cool packs. The product should be stored at 4 °C for immediate use and below -18 °C for long term storage. The product needs to be protected from light.

Instructions for Use

The easy applicable 100x aqueous solution is sterile filtrated and ready-to-go. 1 ml is sufficient for 100 ml of culture medium. One 50 ml bottle will treat 5 liter of culture medium.

ZellShield™ contains a fluoroquinolone and can be used in combination with other antibiotics required as marker. Other components of the product are not applied as genetic markers. The additional use of Pen/Strep is not required or recommended but generally possible. ZellShield™ can also be added directly to any size animal cell culture. Remove and replace fresh ZellShield™ with every time the culture medium is changed. ZellShield™ can be added to freshly passaged and trypsinated cells or to already precultivated cells by medium change or direct application. The reagent should be refreshed by medium change after 7 days of incubation at 37 °C to achieve maximum protection.

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