Mynox[®] Gold

Mycoplasma Elimination Reagent

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
1.1 Kit Components

Instruction Manual

Cat. No. 10-0201 contains 2 packs with Mynox®Gold Reagent

Cat. No. 10-0501 contains 5 packs with Mynox®Gold Reagent

Cat. No. 10-1001 contains 10 packs with Mynox®Gold Reagent

Each pack contains 1 vial *Starter Treatment* and 3 vials *Main Treatment*. Each component is a sterile, ready-to-use solution, aliquoted per vial for single applications of 520 µl/vial.

1.2 Stability and Storage

Mynox®Gold is stable until the expiry date given in the bag label if stored at +2 °C to +8 °C in the dark. Shipping at room temperature.

1.3 Supplemental Requirements

- Standard cell culture equipment and consumables
- Cell culture medium, fetal calf serum, trypsin
- Mycoplasma detection system to verify the elimination success, e.g. Minerva Biolabs's Venor®GeM Mycoplasma PCR Detection Kit

2. Application and Test Principle

Contamination of cell cultures by mycoplasma occurs frequently. For both biological and economical reasons, it is important to eliminate mycoplasma from cell cultures being used for basic research, diagnostics, and biotechnological production. The most commonly used method for elimination, inactivation, or suppression of mycoplasma in cell cultures is treatment with antibiotics. In general, antibiotic therapies alone do not result in long-lasting, successful elimination. Also, the cytotoxic properties of antibiotics can cause undesirable side effects on eukaryotic cells and can facilitate the development of resistant mycoplasma strains.

Mynox®Gold represents the further development of the classical Mynox®. Mynox®Gold is a combination of a standard antibiotic and a biological reagent. In comparison to mammalian cells, mycoplasma lack a cell wall but are encircled by a cytoplasmic membrane. This biological reagent integrates into the mycoplasma membrane and compromises its integrity. By the combination with a standard antibiotic, the effective dose of both, the reagent and the antibiotic, can be reduced to a minimum for lowest cytotoxicity, still causing a highly reliable and definite elimination of mycoplasma. These biophysical properties making the development of resistant strains highly unlikely.

One application comprises of 4 vials, a *Starter Treatment* and three *Main Treatments*. The *Starter Treatment* kills most of the mycoplasma particles without harming the cells. The *Main Treatment* kills all remaining particles leading to a permanent eradication.

Mynox®Gold is intended for research use only.

3. Sample Material

The product is used for the elimination of all kinds of *Mollicutes* and related organisms (e.g. *Mycoplasma*, *Acholeplasma*, *Spiroplasma*, and *Entomoplasma*) in all kinds of cell and virus cultures. For best results the following recommendations should be considered:

1. Not more than 10^5 cells shall be used for a treatment to keep the mycoplasma load low.
2. The mycoplasmacidal activity of Mynox®Gold is affected by the concentration of lipids and proteins in the reaction mixture, e.g. components in fetal calf serum supplement. These ingredients competitively bind Mynox®Gold and prevent its binding to the mycoplasma membrane. A protocol was designed that requires a supplement of 5 % v/v fetal calf serum. Due to the mitigating effect of serum, it is impossible to design a specific protocol that is applicable to the treatment of biologicals with high protein and lipid concentrations.
3. The type of cell culture medium does not effect the efficiency of the treatment. Antibiotics, especially if required for selection, can stay in the elimination mixture. In rare cases cytotoxicity might be increased by unpredictable interactions of the reagents.
4. Viruses shall be treated in combination with their host cells.
5. Mynox®Gold does not penetrate the cellular membrane. Therefore it cannot eliminate intracellular contamination. However, mycoplasma is an extracellular contaminate. *Mycoplasma penetrans* is the only species described intracellularly. Thus far, *M. penetrans* has not been reported as cell contaminant.
6. Since Mynox®Gold works by biophysical means through association with the mycoplasma membrane, the reagent needs direct contact with the mycoplasma particle in order to be effective. Treatment of cell clusters should be avoided. Mycoplasma are protected in intercellular spaces as well as in pockets and clefts of the cell membrane, which can prevent contact with the drug. We suggest using trypsin to detach the cells from each other and to smoothen cell surfaces.

4. Protocol Guide

These protocols have been designed for typical cell lines requiring standard media. The protocol can be used for adherent and suspension cell lines. Minerva Biolabs does not guarantee that these protocols will work in all laboratory situations. Modification of these protocols may be required per individual case.

1. Mix in a sterile culture flasks or petri dish 4.5 ml standard cell culture medium with 5 % v/v FCS and 500 μ l Mynox®Gold Starter Treatment (vial with orange cap)
2. Transfer 5 ml of 1×10^4 to 1×10^5 single cells in cell culture medium with 5% (v/v) FCS into the mix.

The total volume of the treatment mixture is 10 ml.

Ensure the treatment of single cells (check under microscope!). If necessary increase duration of trypsin treatment or detach the cells from each other by pipetting up and down.

Ensure that Mynox®Gold is already present in the culture medium before adding cells. Add cells directly to the elimination mix to avoid evaporation (insert the pipett tip directly into the mix!).


The final FCS concentration is 5 % (v/v).

3. Grow the cell to 80-90% confluency.
4. Split the cells and passage at the usual rate. Add 500 µl of the Mynox®Gold Main Treatment Reagent (vial with transparent cap) to 9.5 ml of passaged cells in fresh media.
Adjust the FCS concentration. Supplementation to 5 % is not required any longer.
5. Repeat Step 3-4 twice. After the third treatment (and a total of 4 passages starting with the starter treatment) the procedure is finished and the culture is free of mycoplasma.

5. Testing for Mycoplasma

Mynox®Gold-treated cell cultures and virus stocks should be subcultivated for four additional passages without mycoplasma-active antibiotics and then assayed for mycoplasma re-emergence to validate culture purity. For highly sensitive detection of mycoplasma contamination, we recommend Venor®GeM Mycoplasma Detection Kit that uses PCR technology (Cat.-No. 11-1025/1050/1100/1250). PCR detection methods should not be used instantly after treatment. Mynox®Gold is lysing the mycoplasma particles and the mycoplasma DNA is released into the culture medium. This DNA would be detected by PCR giving false-positive results. Medium replacement and extracellular DNAses will reduce the level of free mycoplasma DNA within 1 to 2 passages.

Scheme of the treatment of adherent cell lines

Appendix

Limited Product Warranty

This warranty limits our liability for replacement of this product. No warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided. Minerva Biolabs shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Notice to Purchaser: Limited License

The purchase price of this product includes a limited, non-transferable license under German Patent 195 21 938 or its foreign counterpart, owned by Minerva Biolabs GmbH, to use only this amount of the product to practice the Mynox® technique. No right to perform or offer commercial services of any kind using Mynox® is hereby granted by implication or estoppel. Further information on purchasing licenses may be obtained by contacting Minerva Biolabs GmbH, Landhausring 7, 12683 Berlin, Germany.

Trademarks

Mynox, Venor®GeM , Onar and Mycoplasma Off® are registered trademarks of Minerva Biolabs GmbH, Germany.

Minerva Biolabs' International Distributors

Argentina

Genbiotech SRL
Tel: +541145526873
Fax: +541145526873
Email: info@genbiotech.com.ar
Web: www.genbiotech.com.ar

Australia

Biocene Pty. Ltd.
Tel: +61 2 9966 8166
Fax: +61 2 9966 8300
Email: jenny@biocene.com

Austria

BioProducts
Tel: +43 2268 61 65 11
Fax: +43 2268 61 65 44
Email: info@bioproducts.at
Web: www.bioproducts.at

Belgium, Luxembourg

Lucron Bioproducts bvba
Tel: +32 92 82 05 31
Fax: +32 92 82 05 32
Email: info@lucronbioproducts.com
Web: www.lucronbioproducts.com

Canada + USA

Medicorp Inc.
Tel: +1 514 7331 900
Fax: +1 514 7331 212
Email: mktg@medicorp.com
Web: www.medicorp.com

China

Vian-Saga Biological Technology Ltd.
Tel: +86 10 8411 8493
Fax: +86 10 8411 8494
Email: Michael@vian-saga.com

Croatia

Gorea Plus d.o.o.
Tel: +385 1 3369 610
Fax: +385 1 3369 611
Email: management@gorea-plus.hr

Czech Republic

BIO-Consult Laboratoriers spol. sro.
Tel: +42 2 4447 1239
Fax: +42 2 4447 1239
Email: info@bioconsult.cz
Web: www.bioconsult.cz

Finland

Immuno Diagnostic Oy
Tel: +358 3 615 370
Fax: +358 3 682 2039
Email: info@immunodiagnostic.fi
Web: www.immunodiagnostic.fi

France

Biovalley
Tel: +33 1 6007 2020
Fax: +33 1 6007 5051
Email: biovalley@biovalley.fr
Web: www.biovalley.fr

Greece

Bioanalytica S.A.
Tel: +30 210 6400 318
Fax: +30 210 6462 748
Email: bioanalyt@hol.gr
Web: www.bioanalytica.gr

Greece

Varelas S.A.
Tel: +30 210 5281 901
Fax: +30 210 5220 926
Email: varelas@otenet.gr

Great Britain

Cambio Ltd.
Tel: +44 1954 210 200
Fax: +44 1954 210 300
Email: support@cambio.co.uk
Web: www.cambio.co.uk

India

Zelle Biotechnology Pvt. Ltd.
Tel: +91 22 2685 87 41
Fax: +91 22 2685 87 44
Email: anurag@zellebiotech.com

Ireland

Medical Supply Company
Tel: +353 1 8224 222
Fax: +323 1 8224 100
Email: dmcglade@medical-supply.ie
Web: www.medical-supply.ie

Israel

Origolab Ltd.
Tel: +972 2 566 9285
Fax: +972 2 561 2120
Email: origolab@netmedia.net.il

Italy

Biospa
Tel: +39 2 891 391
Fax: +39 2 891 20996
Email: biospa@spaspa.it
Web: www.spaspa.it/biospa.htm

Japan

Funakoshi Co. Ltd.
Tel: +81 3 5684 1615
Fax: +81 3 5684 1775
Email: info@funakoshi.co.jp
Web: www.funakoshi.co.jp

Korea

Bio and Information Corp
Tel: +82 31 7139 439
Fax: +82 31 7139 438
Email: sales@bioninfo.com
Web: www.bioninfo.com

Lithuania

Interlux
Tel: +370 5 27 86 850
Fax: +370 5 27 96 728
Email: marius@interlux.lt
Web: www.interlux.lt

Malaysia

Helix Biotech (M) Sdn Bhd
Tel: +60 390 76 80 10
Fax: +60 390 76 80 07
Email: helixbio@tm.net.my
Web: www.helixbiotech.com

Netherlands

Lucron Bioproducts B.V.
Tel: +31 485 51 16 75
Fax: +31 485 51 20 52
Email: lucron@lucron.nl
Web: www.lucron.nl

New Zealand

Medical & Scientific Ltd.
Tel: +64 96 34 10 36
Fax: +64 96 34 51 46
Email: nzms@nzms.co.nz
Web: www.nzms.co.nz

Norway

E. Pedersen & Son
Tel: +47 22 95 59 59
Fax: +47 22 95 59 40
Email: eped@eped.com
Web: www.eped.com

Poland

STI
Tel: +48 61 641 77 59
Fax: +48 61 641 77 58
Email: office@sti.biz.pl
Web: www.sti.biz.pl

Portugal

Quilaban Lda.
Tel: +21 923 63 50
Fax: +21 923 63 89
Email: quilaban@quilaban.pt
Web: www.quilaban.pt

Slovenia

Kemomed d.o.o.
Tel: +386 4 201 50 50
Fax: +386 4 201 50 55
Email: info@kemomed.si
Web: www.kemomed.si

Spain

LacClinics S.A.
Tel: +34 3 446 47 00
Fax: +34 3 348 10 39
Email: info@labclinics.com
Web: www.labclinics.com

Sweden

ANL-Produkter AB
Tel: +46 8 99 00 90
Fax: +46 8 99 20 40
Email: info@anl.se
Web: www.anl.se

Switzerland

Socochim SA
Tel: +41 21 721 04 50
Fax: +41 21 721 04 51
Email: info@socochim.ch
Web: www.socochim.ch

Taiwan

Only Science Co. Ltd.
Tel: +886 2 2758 5926
Fax: +886 2 2758 6305
Email: cychen@onlyscience.com.tw

Thailand

Biomed Diagnostic Co. Ltd.
Tel: +66 2 8796 026
Fax: +66 2 8796 065
Email: somboon@biomedthai.com
Web: www.biomedthai.com

Turkey

Genomed Saglik Hizmetleri A.S.
Tel: +90 212 248 20 00
Fax: +90 212 220 15 64
Email: muratyazici@genomedtr.com
Web: www.genomed-biotech.com

Related Products for Cell Culture Contamination Control

Diagnostic Kits for conventional PCR

11-1025	Venor®GeM Mycoplasma Detection Kit	25	tests
11-1050	Venor®GeM Mycoplasma Detection Kit	50	tests
11-1100	Venor®GeM Mycoplasma Detection Kit	100	tests
11-1250	Venor®GeM Mycoplasma Detection Kit	250	tests

Diagnostic Kits for real-time PCR

11-4025	Venor®GeM-qEP Mycoplasma ssp. Detection Kit	25	tests
11-4100	Venor®GeM-qEP Mycoplasma ssp. Detection Kit	100	tests
11-4250	Venor®GeM-qEP Mycoplasma ssp. Detection Kit	250	tests
11-6025	Venor®GeM-qDual Mycoplasma Detection Kit with external Inhibition control	25	tests
11-6100	Venor®GeM-qDual Mycoplasma Detection Kit with external Inhibition control	100	tests

Polymerases

53-0050	MB TAQ DNA Polymerase	50	units
53-0100	MB TAQ DNA Polymerase	100	units
53-0200	MB TAQ DNA Polymerase	200	units
53-0250	MB TAQ DNA Polymerase	250	units

Genomic DNA Extracts 100 µl each

51-0101	<i>Legionella pneumophila</i> , ATCC 33152	+/- 10 ng / 100 µl
51-1514	<i>Legionella pneumophila</i> subsp. <i>fraseri</i> , DSMZ 7514	+/- 10 ng / 100 µl
51-1515	<i>Legionella pneumophila</i> subsp. <i>pascuella</i> , DSMZ 7515	+/- 10 ng / 100 µl
51-1368	<i>Legionella bozemanii</i> , NC 011368	+/- 10 ng / 100 µl
51-1370	<i>Legionella dumoffii</i> , NC 011370	+/- 10 ng / 100 µl
51-1723	<i>Lactobacillus acidophilus</i> , NC 001723	+/- 10 ng / 100 µl
51-0566	<i>Streptococcus pneumoniae</i> , DSMZ 20566	+/- 10 ng / 100 µl
51-0792	<i>Clostridium acetobutylicum</i> , DSMZ 792	+/- 10 ng / 100 µl
51-5571	<i>Bordetella pertussis</i> , DSMZ 5571	+/- 10 ng / 100 µl
51-3415	<i>Bordetella parapertussis</i> , DSMZ 13415	+/- 10 ng / 100 µl
51-0177	<i>Ureaplasma urealyticum</i> , NC 10177	+/- 10 ng / 100 µl
51-0111	<i>M. hominis</i> , NC 010111	+/- 10 ng / 100 µl
51-0112	<i>M. orale</i> , NC 010112	+/- 10 ng / 100 µl
51-0113	<i>M. salivarium</i> , NC 010113	+/- 10 ng / 100 µl
51-0115	<i>M. gallisepticum</i> , NC 010115	+/- 10 ng / 100 µl
51-0116	<i>Acholeplasma laidlawii</i> , NC 010116	+/- 10 ng / 100 µl
51-0117	<i>M. fermentans</i> PG19, NC 010117	+/- 10 ng / 100 µl
51-0119	<i>M. pneumoniae</i> , NC 010119	+/- 10 ng / 100 µl
51-0129	<i>M. arginini</i> , NC 010129	+/- 10 ng / 100 µl
51-0130	<i>M. hyorhinis</i> , NC 010130	+/- 10 ng / 100 µl
51-0162	<i>M. arthritis</i> , NC 010162	+/- 10 ng / 100 µl
51-0195	<i>M. genitalium</i> , NC 010195	+/- 10 ng / 100 µl
51-1746	<i>M. penetrans</i> , NC 11746	+/- 10 ng / 100 µl

Quantification Standards 100 µl each

52-0101	<i>Legionella pneumophila</i> DNA Standard	1x10 ⁶	genomes/µl
52-0119	<i>Mycoplasma pneumoniae</i> DNA Standard	1x10 ⁶	genomes/µl
52-0112	<i>Mycoplasma orale</i> DNA Standard	1x10 ⁶	genomes/µl
52-0116	<i>Acholeplasma laidlawii</i> DNA Standard	1x10 ⁶	genomes/µl

Mycoplasma Off®

15-1000	Surface Disinfectant Spray	1000	ml
15-5000	Surface Disinfectant Spray, refill	5000	ml

DNA Remover™

15-2025	DNA Decontamination Reagent, spray bottle	250	ml
15-2200	DNA Decontamination Reagent, refill bottle	4x 500	ml